



Fluorescence cell counter with the best eye for detail

CytoSMART™ Exact FL

Accelerate your research with fluorescence cell counting

Automated fluorescence cell counting has transformed the process of cell quantification, as it offers a great deal of flexibility and accuracy, and can be used when bright-field image analysis is not possible. The CytoSMART Exact FL is an automated, dual fluorescence cell counter with an expanded field of view and unmatched resolution and magnification power. Using an advanced optical system and AI-powered image analysis software, the CytoSMART Exact FL calculates the exact number of cells in a sample and provides reliable assessment of key cellular parameters, including cell viability and transfection efficiency.

The CytoSMART™ Exact FL features:

- + **State-of-the-art optics** allows to accurately distinguish individual cells in clumps and count cells down to 4 µm in diameter
- + **Dual-channel fluorescence** can be used for examining health and viability of hard-to-detect-cells
- + **AI-powered software** performs cell counting, minimizing user-to-user variability
- + **Large field of view & multi-count** let users analyze a larger portion of the available counting volume
- + **Add-ons (e.g. organoid counting)** expand the application range of the device
- + **Reusable or disposable counting slides** for cost- or time-efficient cell counting

Applications

The use of fluorescence cell counting significantly accelerates experimental work in the fields of cancer biology, immunology, tissue engineering, and cell therapy. Numerous cell-based assays significantly benefit from fluorescence cell counting due to the number of issues associated with detecting and analyzing mammalian cells using colorimetric dyes (*e.g.* Trypan Blue). The CytoSMART Exact FL, fitted with red and green fluorescent channels, in combination with appropriate fluorescent dyes, can assist you in examining the following features (Fig. 1):

- + Cell viability
- + Small and hard-to-detect cells
- + Primary cell samples
- + Heterogeneous cell populations
- + Transfection efficiency
- + Distinguishing whole cells from debris
- + and more

State-of-the-art optics

Using the 6.4 MP CMOS camera combined with 10× magnification, the CytoSMART Exact FL can visualize and count cells down to 4 μm in diameter. With increased field of view and multi-count feature, it is also possible to take multiple snapshots of the same sample at various slide positions, which increases the accuracy and reproducibility of a cell count.




Figure 2. The advanced optical system of the CytoSMART Exact FL allows to visualize a range of cells (4 – 70 μm).

AI-powered algorithms

Manual cell counting introduces user-to-user variability to a scientific study. In addition, the process of manual cell counting is challenging and time-consuming, especially when working with cells that are prone to clumping. Equipped with a robust pattern recognition software and a declustering algorithm, the CytoSMART Exact FL detects resuspended single cells, as well as individual cells within cellular clusters with high accuracy and without any bias.




Figure 3. Raw data (left) is processed by the highly trained deep-learning neural network algorithm (right) that can analyze even the most challenging samples.

Automated software updates & secure cloud data storage

All CytoSMART devices are connected to the CytoSMART cloud environment. This means that the data generated by the Exact FL fluorescence cell counter is processed, analysed and stored in one accessible and secure place, protected by a multi-layered cyber security system. Novel software updates are also introduced via the CytoSMART Cloud. With only a few clicks our customers can download and install these add-ons to make their user experience even more seamless. Importantly, the added features will not affect consistency and reproducibility of acquired results.

Innovative organoid counting software

Organoids are stem cell-derived 3D culture systems that have been used to study genetic disorders, infectious diseases, and various types of cancer. To ensure consistency and reproducibility across multiple organoid-based experiments, it is critical to use a fixed concentration of organoids of similar size. Using the Organoid Counting add-on, the CytoSMART Exact FL can analyse bright-field images to determine precise number and size of organoids in a sample.

Frequently Asked Questions

Q: Which fluorescent dyes can be used with the CytoSMART Exact FL?

A: Many different fluorescent dyes can be used, as long as the fluorescent dye's excitation and emission spectra correspond with the fluorescent filters of the Exact FL (green – excitation: 452/45 nm, emission: 512/23 nm; red – excitation: 561/14 nm, emission: 630/90 nm). Some examples are propidium iodide (PI) and red fluorescent protein (RFP) for the red channel, and acridine orange (AO), calcein-AM, and green fluorescent protein (GFP) for the green channel.

Q: Can the CytoSMART Exact FL access cell viability?

A: Yes, the viability of cells can be assessed either using the green/red fluorescent channels (AO/PI) or a bright-field channel (Trypan Blue).

Q: Is a computer required?

A: The device requires a dedicated desktop or laptop running on Windows 10 or above and with a USB 3 port and an active internet connection. A WiFi or Ethernet connection is needed for connecting to the CytoSMART Cloud for image analysis and storage.

Q: Do I need to calibrate the CytoSMART Exact FL?

A: The calibration is not required.

Q: Is it required to use a light cover when performing fluorescence cell counting?

A: To get the best results, we advise covering the stage of the device with the provided stage cover to minimize the amount of environmental light.

Technical Specifications

| | |
|----------------------------|---|
| Counting range | $(5.0 \times 10^4) - (1.0 \times 10^7)$ |
| Size range - cells | 4 – 70 μm |
| Size range - organoids | 10 μm to 100/200 μm^* |
| Loading volume (1 chamber) | 10 μL |
| Field of view | 1 count: 2.1 mm^2 8 counts (multi-count): 16.8 mm^2 |
| Autofocus | Yes |
| Channels | Bright-field, green and red fluorescence channels |
| Fluorescence filters | Green – excitation: 452/45 nm; emission: 512/23 nm Red – excitation: 561/14 nm; emission: 630/90 nm |
| Magnification | 10 \times (digital zoom 20 \times) |
| Image resolution | 2072 \times 2072 pixels |
| Camera | 6.4 MP CMOS |
| Data formats | JPG, CSV, PDF |
| Data storage | The CytoSMART Cloud |
| Computer requirements | Windows 10, USB 3.0 |
| Counting chamber | Reusable (glass) or disposable (plastic) slides |
| Warranty | 1-year and 3-year** |
| Extra features | Live-gating, (pre)dilution calculator, organoid counting * Size limits depend on dimensions of a counting chamber; user can set limits in application ** If a purchase order is made before July 31, 2021 |



(주) 바이오디 www.bio-d.co.kr

(우.14322) 경기도 광명시 하안로 60 광명SK테크노파크 A동 1102호
T_02.6264.3399 F_02.6264.3400 제품문의_cs@bio-d.co.kr 업체문의_ds@bio-d.co.kr

