



XinPOWER-600 Classic Electrophoresis Power Supply

User Manual



Contents

01 Description-----	1
02 Features-----	1
03 Safety Information-----	2
04 Package Contents-----	2
05 Instructions-----	3
06 Specifications-----	5
07 Quality guarantee-----	6
08 Technical Support-----	6

01 Description


Electrophoresis Power Supply is an instrument that helps the electrophoresis cell to realize the normal operation, it mainly consists of power supply, electrophoresis tank, testing unit etc.

XinPOWER-600 adopts the technology of microcomputer and switch, which is more advantageous than the traditional electrophoresis power supply. The main features are time-fixed output and non-time fixed output, stable voltage output, short circuit protection, load-oriented voltage output. It is suitable for the diving-style gel electrophoresis and mini vertical gel electrophoresis.

02 Features

- Constant voltage output, constant current out, constant power output
- Time-set function, automatic switch between constant voltage and constant current in case of exceeding set value.
- Automatic non-load detecting, automatic protection of short circuit, automatic recording.
- 192*64 LCD extra-big screen
- Skid-resistance leg of box

03 Safety Information

 Please read the safety information carefully before use to avoid possible failures which may cause electric shock, and /or damage!

1. Do not touch any part of outlet of electrophoresis power supply during the operation.
 2. Power off the electrophoresis power supply in case of adjusting the cell.
 3. Do not take out the wire or insert the power cord during the operation in order to avoid the damage of the components and protection output.
 4. The power supply design is safe. Please do not disassemble and refit by yourself to avoid safety issues.
 5. We are irresponsible for any harm to the customer caused by any improper repair.
- We are here to inform the customer that there is high voltage inside the main body.

04 Package Contents

- | | |
|--------------------------------|---|
| ● Electrophoresis Power Supply | 1 |
| ● Power Core | 1 |
| ● User Manual | 1 |
| ● Certificate of Approval | 1 |

05 Instructions

1. Key board

There are five buttons equipped with XinPOWER-600 electrophoresis power supply, namely two toggle keys (◀ ▶), two data-adjusting keys (▲ ▼) and one order-implemented key (OK) .

Toggle key: to switch the required voltage.

Data-adjusting key: to adjust the selected data.

Order-implemented key: to implement the already-set order.

Display mode: synchronous voltage, current and power.

Time set: 0-99 hours 59 minutes.

Voltage: 10-600V

Current: 5-1200 mA

Power: 5-500 W

2. Display

XinPOWER-600 Electrophoresis power supply displays the voltage, current, power simultaneously for the real-time checking. The following indicator appears in case of any default:

“NO LOAD” indicates no connecting the load.

“SHORT ” indicates short circuit.

3. The mode of voltage and current can be adjusted

1) Constant voltage output

A. When the icon of electrophoresis power supply appears on the main screen, press the key (▶) to move the cursor to constant value V, at this time, press the key (▲▼) to switch the power supply constant value.

B. While the constant value V appears on the screen, continue to press the key (▶), move the cursor to the specific value of voltage in the middle part of screen, at this time, press the key (▲▼) to adjust the specific value of constant voltage.

C. Continue to press the key (▶) to move the cursor to the position showing the time 00:00, at this time, press the key (▲▼) to set the time.

2) Output of constant current and constant power

3) Operation ends

The power supply to stop output when the operation is finished, then the it is buzzing, press the (▲) key for silent, press OK key to continue operation according to previous set value.

4) Output the record

Press the key (▼) on the screen to enter the page of record, press the key (▲▼) to scroll the record. Press the key of OK to operate according to the value set on the current screen.

06 Specifications

Apparent Parameter	
Unit dimensions	300 mm (L) x 245 mm (W) x 100 mm (H)
Product parameters	
Power supply	100-240V AC, 50/60Hz
Fuse	3A / 250VAC
Maximum power	700W
Output	<ul style="list-style-type: none">• Four output sockets with 3.5mm diameter• Output voltage ripple is less than 1%• Output voltage drift is less than 1% (10 minutes after starting the main body)• Output voltage regulation factor is less than 1%
Output set	<ul style="list-style-type: none">• Output time set range: 0-99 hours 59 minutes, 1-minute adjusting unit• Output constant voltage range: 10-600V, 1 V adjusting unit• Output constant current range: 5-1200mA, 1- mA adjusting unit• Output constant power range: 5-500W, 1W adjusting unit
Safety protection measures	<ul style="list-style-type: none">• Non-load detecting• Protection of short circuit and over-load.
Weight	2.5 Kg
Working conditions	
temperature	0°C-40°C
humidity	0%-90%

07 Quality Guarantee

- (1) The warranty is 1 years since the date of sales.
- (2) The warranty excludes the following situations otherwise it is charged
 - a. Improper operation or accident factors.
 - b. The damage is caused by the user's repair.

08 Technical Support

If you have any questions about the use or operation of the products, please contact us.

Clinx Science Instruments Co., Ltd.

Address: 5C-102, 258 West SongXing Road, Baoshan District, Shanghai, 200940,
China

Email: info@clinxsci.com

Website: www.clinxsci.com