

Shaker for Microplates

MS100



Description:

Features:

- ◆ MS100 Shaker for Microplates is a type of high speed shaker, which is with brushless DC motor and microprocessor controlled technique, specially for elisa plates(96/384 wells), tissue culture plates(24/ 48/96 wells).
- ◆ LCD display system status and parameters.
- ◆ Stable and reliable shaking.
- ◆ Easy to operate with one touch knob.
- ◆ Automatic fault detection and buzzer alarm function.
- ◆ 12-month warranty.

Parameters:

Project	Parameters
Timing range	1min ~ 99h59min
Mixing orbit	3mm
Sample capacity	4plates
Voltage	AC100V~240V, 50/60Hz
Weight(kg)	4.0

Project	Parameters
Shaking speed	200~1500rpm
Speed step value	10rpm
Power	45W
Dimension(mm)	280×270×110
---	---