

Magic Programmable Mixer

360° vertical rotating- for larger particles longer wavelength mixing

Swinging -for culture medium and larger particles medium- longer wavelength mixing

Oscillating- for small particles short wavelength mixing

Feature

- Tube holder assemble discretionary
- Large LCD show the all the parameter
- 2 additional custom programs (6 segments /program) can be created easily
- 10 different programs for rotation and mixing
- 360° rotation and 0-360° discretionary mixing angle is available
- Adjustable speed and time for effective mixing
- Accept tubes from 0.2ml to 50ml
- Cold room and incubator compatible


Possess 10 national patents, use the special antibacterial Materials



The Magic mixer is a combination of variable speed rotator, advanced multiple tube vortex in one machine , it accept tubes from 0.2ml to 50ml ,all kinds of tube holder can be assembled discretionary in one rack.

The mixer can be used in ambient temperatures from +4 to 65°C , making it ideal for a variety of laboratory settings such as Molecular Biology ,Biochemistry ,Clinical and Histochemistry .

Small enough to fit on a crowded bench .the mixer is easily moved where it is needed the bench . cold room ,incubator or biological hood .



Specification

Model	TM-MAGIC MIXER
Adaptor	DC 12V, 1000mA
Rotation speed	1-99rpm
Environment temperature	±4-65°C
Relative air humidity (at 20°C) %	up to 80

Orde information

Cat.No.	Code No	Specification
8011252	TMM-5L	420mm length Magic Mixer with standard rack (D10-12.5mm tubeX20, D16-16.5mm tubeX8,D29-29.5mm tubeX2)
8011252-1	LH5X36	Rack for diameter 10-12.5mm tubes X 36
8011252-2	LH15X22	Rack for diameter 16-16.5mm tubes X 22
8011252-3	LH50X12	Rack for diameter 29-29.5mm tubes X 12
8011252-4	LR15X12	Rack for diameter 16-16.5mm tubes X 12
8011252-5	LR50X12	Rack for diameter 29-29.5mm tubes X 12
8011253	TMM-5M	344mm length Magic Mixer with standard rack (D10-12.5mm tubeX12, D16-16.5mm tubeX2,D29-29.5mm tubeX2)
8011253-1	MH5X24	Rack for diameter 10-12.5mm tubes X 24
8011253-2	MH15X14	Rack for diameter 16-16.5mm tubes X 14
8011253-3	MH50X8	Rack for diameter 29-29.5mm tubesx8
8011253-4	MR15X6	Rack for diameter 16-16.5mm tubes X 6
8011253-5	MR50X6	Rack for diameter 16-16.5mm tubes X 6
8011254	TMM-5S	316mm length Magic Mixer with standard rack (D10-12.5mm tubeX8, D16-16.5mm tubeX2,D29-29.5mm tubeX2)
8011254-1	SH5X20	Rack for diameter 10-12.5mm tubes X 20
8011254-2	SH15X10	Rack for diameter 16-16.5mm tubes X10
8011254-3	SH50X6	Rack for diameter 29-29.5mm tubes X 6
8011254-4	SR15X6	Rack for diameter 16-16.5mm tubes X 6
8011254-5	SR50X6	Rack for diameter 29-29.5mm tubes X 6
8011255	TMM-PL	Mixer 420x125x168 with Platform (267X154)
8011256	TMM-PM	Mixer 344x125x168 with Platform (169X154)
8011257	TMM-PS	Mixer 316x125x168 with Platform (139X154)