

Mivot mixer

The Vortex mixer widely used in life science, physical and chemical analysis fields, for vigorous re-suspension of cells or chemical pellets in tubes up to 15 ml.

Features

Touch and continuous operations

Mini and compact appearance , powerful function , sucker type under-chassis , super shockproof, ideal for high speed operation.

Wide speed range:0-3000rpm/min, continuously variable motor and Brushless DC Motor are optional. Stable low speed, powerful high speed.

Shaking block convenient for replacement, stable and reliable, sturdy and durable eccentric bearing design.

Various shaking block adapters optional.(Suitable for Eppendorf Tube)



Specification

Orbit	4mm
Shaking mode	circle
Operation mode	Continuous running/ Short Mix
Speed Range	100 ~ 3000rpm
Voltage	AC220V/AC110V,50/60Hz
Power	60W



MT-A



MT-B



MT-C



MT-D



MT-E



MT-F



MT-G

Order information

Cat.No.	Part.No.	Specification
8011690	TM-Mivot	Vortex Mixer 220V/50Hz(with standard head)
8011690-1	MT-S	Standard head. Suit for 1.5mlX8 centrifuge tubes
8011690-2	MT-A	Foam tube holder for 15mlX8 & 10mlX8
8011690-3	MT-B	50mlX6/Foam tube holder for 50mlX6
8011690-4	MT-C	Plate holder platform. For all standard plates
8011690-5	MT-D	Plate adapter suitable for less than 120mm diameter tube or container .
8011690-6	MT-E	Tube holder platform. For 1.5ml/2.0ml tubes x 12
8011690-7	MT-F	Tube holder platform. For 15ml tubes x 4
8011690-8	MT-G	Tube holder platform. For 50ml tubes x 2