

Specifications	XinPOWER-300 Basic Power Supply	XinPOWER-600 Universal Power Supply
Type of output	constant voltage,constant current,constant power with automatic crossover	constant voltage,constant current,constant power with automatic crossover
Output range (programmable)	0-99 hours 59 minutes: adjustable in one-minute steps 10-300V: adjustable in 1V steps 5-400mA: adjustable in 1mA steps 5-75W: adjustable in 1W steps Tip: The current may be greater than 400mA when the four sets outputting at the same time.	0-99 hours 59 minutes: adjustable in one-minute steps 10-600V: adjustable in 1V steps 5-1200mA: adjustable in 1mA steps 5-500W: adjustable in 1W steps
Programmable methods	Stores 20 methods,each with up to 10 steps	Stores 20 methods,each with up to 10 steps
Display	192*64 Pixel with LCD	192*64 Pixel with LCD
Timer control	1min-99 hr 59 min,fully adjustable	1min-99 hr 59 min,fully adjustable
Safety features	No-load detection,overload/short circuit protection,sudden load change detection,auto power-up after power failure	No-load detection,overload/short circuit protection,sudden load change detection,auto power-up after power failure
Safety fuse	3A/250VAC	3A/250VAC
Power supply	100-240VAC,50/60Hz	100-240VAC,50/60Hz
Maximum power	100W	700W
Operating conditions	0°C-40°C,0%-90% humidity	0°C-40°C,0%-90% humidity
Dimension (L×W×H)	300 mm x 245 mm x 100 mm	300 mm x 245 mm x 100 mm
Net Weight	2kg	2.5Kg
Product picture		