

Bead Ruptor 12

BEAD MILL HOMOGENIZER



*Compelling performance
at a compelling price*

 **OMNI**
International
The Homogenizer Company®

Bead Ruptor 12



ORDER NO. 19-050A

Ideally suited for nucleic acid, drug, protein, enzyme, and virus extractions.

Specifications

Processing Range: 25 μ L to 6 mL

Power Rating: 115-230 VAC/50-60 Hz

Dimensions: 13" W x 15.5" D x 13" H

Weight: 29 kg (65 lbs.)

Compliance: CE

Run Time: 1s - 9:59 m

Performance Range: 0.8 m/s - 6 m/s
in increments of 0.15 m/s

Cycles: 1-10

Includes Tube Carriage for 2 mL and 7 mL tubes

- Holds 12 x 2 mL or 0.5 mL tubes and/or 4 x 7 mL tubes
- Includes semi-permanent ring



Features

- Optimized tube motion with speeds up to 6 m/s
- Process up to 16 samples per cycle
- Convenient front loading design
- Locking safety lid
- Reproducible extraction of proteins, nucleic acids, and small molecules

