

Bead Ruptor 4

BEAD MILL HOMOGENIZER



*Powerful performance
in a compact design*



OMNI
International
The Homogenizer Company®

Bead Ruptor 4

Built upon the innovative design and unsurpassed quality of Omni International's Bead Ruptor family.

Specifications

Processing Range: 25 μ L to 5 mL

Power Rating: 100-230 VAC/50-60 Hz

Dimensions: 8.5" W x 11.5" D x 10" H

Weight: 16.5 lbs. (7.48 kg)

Run Time: 1 - 5 mins.

Programs: 1

Cycles: 1

Performance Range: 1 - 5 m/s

Omni Bead Ruptor 4 is ideally suited for:

- Homogenizing up to 4 samples in seconds
- Reproducible extraction of proteins and nucleic acids
- Simultaneously processing 4 x 0.5 mL, 4 x 1.5 mL, 4 x 2 mL or 1 x 7 mL (sold separately)

Efficiently lyse:

- Human tissues
- Animal tissues
- Plant & food
- Micro-organisms




ORDER NO. 25-010

Designed for personal use

- Broad performance range — gentle mixing to high force homogenization
- Small footprint — requires minimal bench space
- Holds 4 x 0.5 mL, 1.5 mL, 2 mL tubes or 1 x 7 mL (sold separately)

Patent Pending Snap In Tube Carriage



For best results, use with Omni International pre-filled bead kits



.5 mL

1.5 mL

2 mL