

# Tissue Master

## TISSUE HOMOGENIZER

The **Tissue Master 125** is an extremely economical homogenizer that is well suited for customers that are occasionally homogenizing small volume samples with little variation in sample size. The Tissue Master 125 comes complete with a 7 mm stainless steel probe.

The **Omni Tissue Master 125 Homogenizer Complete Kit** includes a 5 mm, 7 mm and 10 mm stainless steel generator probes, specifically designed for processing volumes of <200  $\mu$ L to 100 mL. These probes can handle a wide variety of different sample sizes, making the **Tissue Master** a great choice for a laboratory with many applications, but a limited budget.

**Processing Range:** <200  $\mu$ L to 100 mL

**Speed:** Variable from 500 to 35,000 rpm

**Power Rating:** 125 watt

**Weight:** 535 g (18.8 oz.)

**Sound Level:** <72 db

**Warranty:** 1 year motor warranty

**Standards Approval:** CE Approved

### Features

- Powerful 125 watt motor
- Light weight ergonomic design
- Value priced, complete with 7 mm probe
- Complete kit includes 5 mm, 7 mm and 10 mm stainless steel generator probes



# Omni Stainless Steel generator probes

## OMNI TISSUE MASTER 125 HOMOGENIZER COMPLETE KIT

(WITH 5, 7 AND 10 MM PROBES) ALSO AVAILABLE INDIVIDUALLY

Available in 115V or 220V

ORDER NO. TM125-115

TM125-220



5 MM PROBE

7 MM PROBE

10 MM PROBE

## VALUE PRICED HOMOGENIZING KIT

WITH 5 MM PROBE    **ORDER NO. TM125-115-05**

WITH 7 MM PROBE    **ORDER NO. TM125-115**

WITH 10 MM PROBE    **ORDER NO. TM125-115-10**